


Neuroimaging: basic workflow, visualisation & clinical tractography

Klinik für Neurochirurgie

Dr. Effy Ntemou
Research group leader
Cognitive Neurosurgery Düsseldorf
effrosyni.ntemou@med.uni-duesseldorf.de



File formats

DICOM vs NIfTI

.dcm


- What the scanner generates
 - Clinical standard
 - 1 file = 2D slice
- NOT compressed/ cannot be compressed
 - Big data file
- Works with PACS and RIS
- Metadata = contain personal info!

.nii/.nii.gz

- Research standard
- Stores entire 3D/4D volumes in ONE file
- Easier for processing and sharing
- Metadata = orientation + voxel size
 - Allows for lossless compression (.nii.gz)
 - NOT compatible with PACS
- But with every other neuroimaging software



File formats



...	Test	>	DICOM	>	0000F778	>	AA60E917	>	AA6586E0	>	AA6586E0 durchsuchen	🔍
📁	📁	🗑️	↕ Sortieren	☰ Anzeigen	...	Details						
Name	Änderungsdatum	Typ	Größe									
📁 0000A74	11.11.2025 14:52	Dateiordner										
📁 0000A106	11.11.2025 14:52	Dateiordner										
📁 0000AEDC	11.11.2025 14:52	Dateiordner										
📁 0000B23B	11.11.2025 14:52	Dateiordner										
📁 0000BEAD	11.11.2025 14:52	Dateiordner										
📁 0000CFFC	11.11.2025 14:52	Dateiordner										
📁 0000E6C3	11.11.2025 14:52	Dateiordner										
📁 0000E966	11.11.2025 14:52	Dateiordner										
📁 0000F27A	11.11.2025 14:52	Dateiordner										
📁 0000F31C	11.11.2025 14:52	Dateiordner										
📁 000002AA	11.11.2025 14:52	Dateiordner										
📁 00005B03	11.11.2025 14:52	Dateiordner										
📁 00007E95	11.11.2025 14:52	Dateiordner										
📁 00008A45	11.11.2025 14:52	Dateiordner										
📁 0000015F	11.11.2025 14:52	Dateiordner										
📁 000026B9	11.11.2025 14:52	Dateiordner										
📁 000045A4	11.11.2025 14:52	Dateiordner										
📁 0000121B	11.11.2025 14:52	Dateiordner										
📁 00001600	11.11.2025 14:52	Dateiordner										

How do I find T1/T2/DWI??

The best way is to look at the **number of slices/files**.

Given our current protocol...

T1 (without KM) = 176 slices

T1 (with KM) = 176 slices

T2 = 30

T2 FLAIR = 35

File formats

0000F27A					000002AA				
Presentationen					Presentationen				
0000F778 > AA60E917 > AA6586E0 > 0000F27A					0000F778 > AA60E917 > AA6586E0 > 000002AA				
Neu Sortieren Anzeigen					Neu Sortieren Anzeigen				
Name	Änderungsdatum	Typ	Größe		Name	Änderungsdatum	Typ	Größe	
EE0A8E2D	11.11.2025 14:52	Datei	201 KB		EE2CA41B	11.11.2025 14:52	Datei	567 KB	
EE0AC5B3	11.11.2025 14:52	Datei	201 KB		EE2CF247	11.11.2025 14:52	Datei	567 KB	
EE0D05E6	11.11.2025 14:52	Datei	201 KB		EE5C9474	11.11.2025 14:52	Datei	567 KB	
EE0E5A62	11.11.2025 14:52	Datei	201 KB		EE7A4504	11.11.2025 14:52	Datei	567 KB	
EE0E440D	11.11.2025 14:52	Datei	201 KB		EE7C2F7C	11.11.2025 14:52	Datei	567 KB	
EE1A5757	11.11.2025 14:52	Datei	201 KB		EE7DB9E3	11.11.2025 14:52	Datei	567 KB	
EE1CD957	11.11.2025 14:52	Datei	201 KB		EE09C265	11.11.2025 14:52	Datei	567 KB	
EE1F3E64	11.11.2025 14:52	Datei	201 KB		EE50EDB0	11.11.2025 14:52	Datei	567 KB	
EE2C4C50	11.11.2025 14:52	Datei	201 KB		EE50FEFD	11.11.2025 14:52	Datei	567 KB	
EE2C8839	11.11.2025 14:52	Datei	201 KB		EE052BC7	11.11.2025 14:52	Datei	567 KB	
EE2DA422	11.11.2025 14:52	Datei	201 KB		EE88A220	11.11.2025 14:52	Datei	567 KB	
EE3BB625	11.11.2025 14:52	Datei	201 KB		EE98B0F6	11.11.2025 14:52	Datei	567 KB	
EE3D881C	11.11.2025 14:52	Datei	201 KB		EE367E90	11.11.2025 14:52	Datei	567 KB	
EE3FC7DC	11.11.2025 14:52	Datei	201 KB		EE755A1D	11.11.2025 14:52	Datei	567 KB	
EE4A090B	11.11.2025 14:52	Datei	201 KB		EE807DBE	11.11.2025 14:52	Datei	567 KB	
EE4DDD01	11.11.2025 14:52	Datei	201 KB		EE2633A4	11.11.2025 14:52	Datei	567 KB	
EE4DFDCE	11.11.2025 14:52	Datei	201 KB		EE73013E	11.11.2025 14:52	Datei	567 KB	
EE4EB236	11.11.2025 14:52	Datei	201 KB		EEA02E42	11.11.2025 14:52	Datei	567 KB	
EE4ECC6B	11.11.2025 14:52	Datei	201 KB		EE4E5A7C	11.11.2025 14:52	Datei	567 KB	
176 Elemente					30 Elemente				

File formats

Converting between .dcm and .nii

There are scripts that you can run:

Python workflow:


```
import dicom2nifti
dicom2nifti.convert_directory('/path/to/dicom_folder', '/output/folder')
```

Bash workflow (comes with FSL, SPM, etc.):

```
dcm2niix -o /output/folder -f pat01_T1 /path/to/dicom_folder
```

GUI version (essentially runs dcm2niix):

MRicroGL (<https://www.nitrc.org/projects/mricrogl/>)




If you write your thesis with me, this program will likely become one of your best friends.

It has everything a “good friend” should have:



- It is free and open-access
- It is user-friendly and comes with a GUI
- It is great for visualization
- It runs on every OS

To convert a .dcm to .nii:

- “Import”
- “Convert DICOM to NIfTI”
- You should get something like this
- Take the folder that contains all DICOMs you want to convert and drag and drop it in this window



If it worked properly, you should see something like this.




You can now open your nifti by going to “File → Open”.

The goal now is to **segment** the lesion.

To do this, we need to “Draw” the lesion.

What does this segmentation really create??



Semi-automatic lesion segmentation with MRlcroGL

- If you segmented the lesion every 2 slices, the sagittal plane should look like this.


- **Beginning** and **end** are very important.

- Now we need to turn this into a mass and that's why we need to interpolate.

- Draw > Advanced > Interpolate between slices > All axial gaps

- Now you shouldn't see any gaps along the axial plane.

- Sometimes, lesion edges are rather "sharp" → smooth




Semi-automatic lesion segmentation with MRlcroGL

To save your wonderful creation:


Draw > Save VOI > save it as .nii or .nii.gz

If we open the lesion file using MRlcroGL (as the main image), what will it look like?




Spaces

Individual vs Normalized space



My brain

Location	
Coordinates: Scanner anatomical	
History Location	-2.573767
	55.75864
	33.5925
Volume	




MNI152

Location	
Coordinates: MNI152	
History Location	-0.742368
	32.35208
	4.126101
Volume	

Normalisation


How can my brain fit an average brain?

Most importantly for us, how can we make a lesion fit an average brain?



Rorden & Brett, 2000; *Behav. Neurol.*



Linear transformation
- Affine (preserves straight lines and proportions)



Non linear transformation
- Based on curvature (gyri/sulci)



Normalisation with SPM25

- You can minimize the graphic and processes windows




Normalisation with SPM25


- Go where you downloaded the dummy data and double-click on “defaced_test1.nii.gz”
- You see a new file created with the name “defaced_test1.nii”



Normalisation with SPM25



- Go to your second batch (i.e., Normalise: Write) on the “Module List”
- Double-click on “Data”




Normalisation with SPM25

Here, choose the lesion mask you created before!

For “Interpolation” choose „Nearest Neighbour“.

Choose a filename prefix that is informative to you. I went with “norm_”.




Normalisation with SPM25

Add another “Normalise: Write” module.

Last one, I promise.

For Images to Write under Data select “defaced_test1.nii”



Choose a file prefix that is informative to you. I went for “norm ”.


If everything is ready, click on the little “Play” icon to run our Batch script!

Be a bit patient...

Normalisation with SPM25


What was the output of the 3 modules we ran?

How can we check whether our normalization was successful?



Normalisation with SPM25

Which structures are affected by the tumor?



Descriptives for "natbrainlab" (Copied to clip board)


Background image: natbrainlab

Regions:

TotalVoxelsNotZero 3010


Background image: natbrainlab

Index	Name	numVox	numVoxNotZero	fracVoxNotZero
2	"Arcuate_Anterior_Segment_Left"	4095	2366	0.577777802944183
3	"Long_Segment_Left"		2984	0.21581768989563




Lesion overlap figures

Let's say you have several normalised lesions and you would like to know where the maximum overlap is. How can you do that?




Brainlab & clinical tractography

Choose a patient and then → Cranial Planning



Brainlab & clinical tractography



Muller et al., 2018; *JNS*


Brainlab & clinical tractography

2 ROI approach: 1 seeding + 1 refining

CST

Seed ROI: Brainstem

When you add the seed ROI you get something like this...




Brainlab & clinical tractography

CST

Refining: Precentral

You can split Leg from Hand
CST.

The articulatory CST is
trickier than leg and hand.




Brainlab & clinical tractography

Arcuate Fasciculus

Create a Fiber Bundle

Seeding ROI:
Create ROI

- On the sagittal plane set your crosshair as medial so that you see the beginning of the STS.
- Then go to the axial plane and create a circle with ca. 1.5 -3cm radius around the center of the crosshair.
- “Add” your seeding ROI to the Fiber Bundle.




Brainlab & clinical tractography

Arcuate Fasciculus

Additional ROI:

- On the sagittal plane, find the deepest point of the post/central sulcus. Scroll to move the crosshair as medial so that you can see the beginning of the STS.
- Create a circular ROI with approx. 1.5 – 3cm radius around the crosshair.



Color: Blue



Brainlab & clinical tractography

Arcuate Fasciculus

Now add your refining ROI to the Fiber Bundle!




If you are not totally happy, you can edit the Fiber Bundle by clicking “Remove Fibers”!

Brainlab & clinical tractography

Inferior Longitudinal Fasciculus (both ROIs on coronal)

Seed: Anterior temporal lobe, coronal plane




Brainlab & clinical tractography

Inferior Longitudinal Fasciculus

Refining: Occipital lobe, posterior right above the “end” of the cerebellum. Coronal plane.


Color: Purple



Brainlab & clinical tractography

Inferior Fronto-Occipital Fasciculus (both ROIs on coronal)

Seed: Prefrontal cortex.




Brainlab & clinical tractography

Inferior Fronto-Occipital Fasciculus (both ROIs on coronal)

Refining: You can reuse the ILF refining ROI. These two tracts terminate in very similar occipital regions.


Color: Green



Brainlab & clinical tractography

Frontal Aslant Tract (Axial + Sagittal)

Seed: Place a circle with a 1.5-3cm radius approx. over the SMA using the coronal plane.




Brainlab & clinical tractography


Frontal Aslant Tract (Axial + Sagittal)

Refining: Broca area ROI using the sagittal plane to capture IFG terminations.

Color: Pink



Brainlab & clinical tractography





For questions/comments:
effrosyni.ntemou@med.uni-duesseldorf.de